

Innovators in Microwave Technology
3100 Smith Farm Road, P.O. Box 200, Matthews, NC 28106-0200 USA
Phone (704) 821-7015 * Fax (704) 821-5185 * Email: sales@cem.com * Web: http://www.cem.com

| QUOTATION | | | |
|--------------|--|--|-----------------------------|
| Name | Dr. Maciej Stawikowski | Created Date | 11/11/2015 |
| Phone | (561) 799-8141 | Quote Number | 00008691 |
| Email | mstawikowski@fau.edu | Expiration Date | 1/8/2016 |
| | | | |
| Bill To Name | Florida Atlantic University | Account Name | Florida Atlantic University |
| Bill To | 5353 Parkside Dr Jupiter, FL 33458 United States | Please Provide "Billing and/or Ship To" Information If Different Than Quote: | |
| | | Contact Address City Zip | State/Province |
| | | Phone | Fax |



Innovators in Microwave Technology
3100 Smith Farm Road, P.O. Box 200, Matthews, NC 28106-0200 USA
Phone (704) 821-7015 * Fax (704) 821-5185 * Email: sales@cem.com * Web: http://www.cem.com

| Product Code | Product | Product Description | Sales Price | Quantity | Total Price |
|--------------------------------|--|--|--------------|----------|--------------|
| 925115 | Discover LabMate™ w/ IntelliVent™ Technology, 120 V | Features: -Ample power (0-300W) for a broad range of syntheses -IntelliVent™ pressure device to manage pressure & provide safe venting if over pressurization occurs -Max. operating pressure: 21bar (300psi) -Runs pressurized rxns in 10mL vial (no tools required to assemble) -Runs reflux rxns using standard lab glassware (5-125mL RBF) -Large cavity, w/ protective spill cup, grants full-hand access to simplify routine maintenance -Volume-independent IR temp control (0-300°C) reduces # of vial sizes to stock -Variable speed electromagnetic stirring -Utilizes compressed air to rapidly cool (quench) rxns -System control from easy-to-use Synergy™ software or on-board keypad -Communication ports for accessories & external control -Standard electrical outlet required Included: 10mL vials/caps (100 each), open vessel kit, stir bars, Microwave Synthesis book (CD-ROM), Synergy™ software, & operations manual Available Options: 80mL pressure vessel (scale up); Gas Addition Kit (rxns w/ gaseous reagents); Explorer24 automation deck; CoolMate™ (sub-ambient chemistries); Upgradeable for peptide synthesis &/or enzymatic digest | \$21,500.00 | 1.00 | \$21,500.00 |
| CEM In-Kind Contribution | CEM In-Kind Contribution - Discover | | (\$3,225.00) | 1.00 | (\$3,225.00) |
| 925605 | Liberty Blue Upgrade | The Liberty Blue™ Upgrade converts a Discover to an automated microwave peptide synthesizer designed for Fmoc solid phase peptide synthesis. The system is capable of automated synthesis of a single peptide. The standard synthesis scale is 0.1 to 0.5 mmol using a 30mL reaction vessel. The Liberty Blue™ uses a Discover single mode microwave reactor that provides for self-tuning with a continuous power delivery from 0-300 Watts in 1-Watt increments. The system contains 20 ports standard for amino acids, 7 additional unusual amino acid ports, and 4 external ports for other reagents. All reagent deliveries are controlled by nitrogen pressure. Amino acid, activator, and activator base are delivered with the Flex-Add technology allowing for precise and accurate delivery. Includes Notebook Computer with Application Software and Installation & Training - P/N 923000BS | \$66,500.00 | 1.00 | \$66,500.00 |
| CEM In-Kind Contribution | CEM In-Kind Contribution - Liberty | | (\$9,975.00) | 1.00 | (\$9,975.00) |

Sales Rep Greg LeBlanc **Total Price** \$74,800.00

Currency & Shipping Terms



Innovators in Microwave Technology

3100 Smith Farm Road, P.O. Box 200, Matthews, NC 28106-0200 USA Phone (704) 821-7015 * Fax (704) 821-5185 * Email: sales@cem.com * Web: http://www.cem.com

Standard Delivery is 28 days ARO

| Purchase Order | | | | | | | |
|--|-------|------------|--|--|--|--|--|
| Note: This quote can be used as a purchase order by entering purchase order number and signing | | | | | | | |
| Please sign here: | Date: | PO Number: | | | | | |
| Tax Status: If Account is Tax-Exempt, please provide a Tax-Exemption or Resale Certificate along with the Purchase Order. Otherwise, the order will be processed as taxable. | | | | | | | |
| Standard Terms & Conditions | | | | | | | |

For information on leasing CEM products, contact Marissa Reinhardt at CEM Financial Services Phone 973-292-0025 x412, mreinhardt@captivelease.com

Standard Terms & Conditions

1. Payment - Make all checks payable to CEM Corporation and mail to Accounts Receivable at the below address:

12750 COLLECTIONS CENTER DRIVE, CHICAGO IL 60693 United States

- 2a. (USA) General Net 30 Days FOB to Matthews, NC. Freight is prepaid and added to invoice. Clerical errors are subject to corrections.

 2b. (Canada) General Net 30 Days; Canada Shipping Terms: ExWorks PrePaid & Charge EXP Matthews, NC USA. If DDP Shipping Terms are required, then contact CEM Corporation for revised quotation. CEM Corporation does not withold nor remit GST or QST taxes for Canada. Customer is responsible for remitting these directly. Clerical errors are subject to corrections.
- 3. All orders subject to shipping and handling charges.
- 4. Taxes All prices are quoted excluded of any sales, excise, use, or similar tax. The amount of tax applicable is the buyers responsibility.
- 5. Renters Casualty Renter shall bear the risk of loss or damage to the equipment from delivery to customer site until it is returned to CEM.
- 6. Warranty All instruments are warrantied against defects in workmanship or material for one (1) year from date of shipment.
- 7. Returns All returns must be authorized with an RMA# by the CEM Customer Service Department. Call (800) 726-3331 for authorization.
- 8. Quantity Any changes to quantities other than complete instruments will not effect the validity of this quote.